FLEX
Monoclonal Mouse
Anti-Human
Desmin
Clone D33
Ready-to-Use
(Dako Autostainer/Autostainer Plus)

Code IS606

ENGLISH	
Intended use	For in vitro diagnostic use. FLEX Monoclonal Mouse Anti-Human Desmin, Clone D33, Ready-to-Use, (Dako Autostainer/Autostainer Plus), is intended for use in immunohistochemistry together with Dako Autostainer/Autostainer Plus instruments. This antibody labels smooth and striated muscle cells as well as mesothelial cells, and is useful for the identification of rhabdomyosarcomas (1, 2), leiomyomas, (3, 4) and mesotheliomas (5). The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Summary and explanation	Desmin belongs to the class III of intermediate filaments, constituting part of the cytoskeleton, and is the characteristic intermediate filament of all three types of muscle cells (skeletal, cardiac and smooth muscle). Desmin is a 53-kDa protein encoded by nine exons of a gene located at 2q35. Desmin forms a cytoskeletal network across the muscle fibre bordering at the plasma and nuclear membrane and is particularly localized to the subplasmalemmal region and the Z-band (6). Refer to Dako's <i>General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.
Reagent provided	Ready-to-use monoclonal mouse antibody provided in liquid form in a buffer containing stabilizing protein and 0.015 mol/L sodium azide. <u>Clone</u> : D33. <u>Isotype</u> : IgG1, kappa.
Immunogen	Desmin purified from human muscle.
Specificity	In Western blotting of cytoskeletal preparations of MCF-7 cells, SCC-4 cells RT-112 cells, skin epidermis, primary culture of renal carcinoma cells, leiomyoma, fetal heart, and crude extract from human white matter of brain, the antibody labels only a ~53 kDa band in the leiomyoma and fetal heart preparations, proving the reactivity with desmin, and demonstrating the absence of reactivity with other types of intermediate filaments tested (7). As demonstrated by immunohistochemistry on formalin-fixed, paraffin embedded tissue sections, the antibody cross-reacts with the desmin-equivalent protein in mouse and rat.
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used. 4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 5. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
Specimen preparation including materials required but not supplied	The antibody can be used for labeling formalin-fixed, paraffin-embedded tissue sections. Tissue specimens should be cut into sections of approximately 4 µm. Pre-treatment with heat-induced epitope retrieval (HIER) is required using Dako PT Link (Code PT100/PT101). For details, please refer to the PT Link User Guide. Optimal results are obtained by pretreating tissues using EnVision™ FLEX Target Retrieval Solution, High pH (50x) (Code K8010/K8004). <u>Paraffin-embedded sections</u> : Pre-treatment of formalin-fixed, paraffin-embedded tissue sections is recommended using the 3-in-1 specimen preparation procedure for Dako PT Link. Follow the pre-treatment procedure outlined in the package insert for EnVision™ FLEX Target Retrieval Solution, High pH (50x) (Code K8010/K8004). Note: After staining the sections must be dehydrated, cleared and mounted using permanent mounting medium. <u>Deparaffinized sections</u> : Pre-treatment of deparaffinized formalin-fixed, paraffin-embedded tissue sections is recommended using Dako PT Link and following the same procedure as described for paraffin-embedded sections. After staining the slides should be mounted using aqueous or permanent mounting medium. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. For greater adherence of tissue sections to glass slides, the use of FLEX IHC Microscope Slides (Code K8020) is recommended.
Staining procedure including materials required but not supplied	The recommended visualization system is EnVision™ FLEX, High pH, (Dako Autostainer/Autostainer Plus) (Code K8010). The staining steps and incubation times are pre-programmed into the software of Dako Autostainer/Autostainer Plus instruments, using the following protocols: Template protocol: FLEXRTU2 (200 µL dispense volume) or FLEXRTU3 (300 µL dispense volume) Autoprogram: Desmin (without counterstaining) or DesminH (with counterstaining) The Auxiliary step should be set to "rinse buffer" in staining runs with ≤10 slides. For staining runs with >10 slides the Auxiliary step should be set to "none". This ascertains comparable wash times. All incubation steps should be performed at room temperature. For details, please refer to the Operator's Manual for the dedicated instrument. If the protocols are not available on the used Dako Autostainer instrument, please contact Dako Technical Services. Optimal conditions may vary depending on specimen and preparation methods, and should be determined by each individual laboratory. If the evaluating pathologist should desire a different staining intensity, a Dako Application Specialist/Technical Service Specialist can be contacted for information on re-programming of the protocol. Verify that the performance of the adjusted protocol is still valid by evaluating that the staining pattern is identical to the staining pattern described in "Performance characteristics". Counterstaining in hematoxylin is recommended using EnVision™ FLEX Hematoxylin, (Dako Autostainer/Autostainer Plus) (Code K8018). Non-
(115224-002)	IS606/EFG/MNI/2009.12.04 p. 1/4

Dako Denmark A/S | Produktionsvej 42 | DK-2600 Glostrup | Denmark | Tel. +45 44 85 95 00 | Fax +45 44 85 95 95 | CVR No. 33 21 13 17

aqueous, permanent mounting medium is recommended.

Positive and negative controls should be run simultaneously using the same protocol as the patient specimens. The positive control tissue should include colon and the cells/structures should display reaction patterns as described for this tissue in "Performance characteristics" in all positive specimens. The recommended negative control reagent is FLEX Negative Control, Mouse, (Dako Autostainer/Autostainer Plus) (Code IS750).

Staining interpretation
Cells labeled by the antibody display cytoplasmic staining. The staining may show a fibrillar aspect (7).

Performance characteristics

Normal tissues: The antibody labels vascular smooth muscle cells and esophageal smooth muscle cells (3). In 8/10 fixed tissues, the antibody labeled mesothelial cells (5). Except for cells of muscle origin, normal bone marrow, kidney, breast, prostate gland, brain (gray matter), jejunum, thyroid gland, lung, colon, stomach, ureter and tonsil are negative (1). The smooth muscle cells in the tunica muscularis in colon show a moderate to strong staining reaction, whereas smooth muscle cells in the small vessels in the lamina propria show a weak to moderate staining intensity.

Abnormal tissues: Of human rhabdomyosarcomas, the antibody labeled 4/4 tumors (1). In the human esophagus the antibody labeled 35/35 leiomyomas, 1/3 leiomyosarcomas and 3/17 stromal tumors (3). In addition, the antibody labeled 33/33 highly cellular leiomyomas of the uterus (4). Also, 8/16 malignant, diffuse mesotheliomas of the pleura, were labeled by the antibody (5). In human vaginal and endocervical polyps, the antibody labeled both background stromal cells and the atypical stromal cells within the subepithelial zone (8). The antibody showed no labeling in 10 cases each of neuroblastoma, lymphoma, Wilms' tumor, primitive neuroectodermal tumor, and retinoblastoma (2), nor did it label 4 cases of small-cell carcinoma of lung, 3 cases of Ewing's sarcoma (1), and 6 cases of endometrial stromal nodules (4).

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1.

Chang TK, Li CY, Smithson WA. Immunocytochemical study of small round cell tumors in routinely processed specimens. Arch Pathol Lab Med 1989; 113:1343-8.

2.

Pollock L, Rampling D, Greenwald SE, Malone M. Desmin expression in rhabdomyosarcoma: influence of the desmin clone and immunohistochemical method. J Clin Pathol 1995; 48:535-8.

3.

Miettinen M, Sarlomo-Rikala M, Sobin LH, Lasota J. Esophageal stromal tumors. A clinicopathologic, immunohistochemical, and molecular genetic study of 17 cases and comparison with esophageal leiomyomas and leiomyosarcomas. Am J Surg Pathol 2000; 24:211-22.

4.

Oliva E, Young RH, Clement PB, Bhan AK, Scully RE. Cellular benign mesenchymal tumors of the uterus. A comparative morphologic and immunohistochemical analysis of 33 highly cellular leiomyomas and six endometrial stromal nodules, two frequently confused tumors. Am J Surg Pathol 1995; 19:757-68.

5.

Hurlimann J. Desmin and neural marker expression in mesothelial cells and mesotheliomas. Hum Pathol 1994; 25:753-7.

6.

Goebel HH, Warlo IAP. Progress in desmin-related myopathies (review). J Child Neurol 2000; 15:565-72.

7.

Van Muijen GNP, Ruiter DJ, Warnaar SO. Coexpression of intermediate filament polypeptides in human fetal and adult tissues. Lab Invest 1987; 57:359-69.

8.

Pitt MA, Roberts ISD, Agbamu DA, Eyden BP. The nature of atypical multinucleated stromal cells: a study of 37 cases from different sites. Histopathol 1993; 23:137-45.

Explanation of symbols/ Explication des symboles/ Erläuterung der Symbole			
REF	Catalogue number Référence catalogue Bestellnummer	2 °C  8 °C Temperature limitation Limites de température Zulässiger Temperaturbereich	 Use by Utiliser avant Verwendbar bis
IVD	In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum	 Contains sufficient for <n> tests Contenu suffisant pour <n> tests Inhalt ausreichend für <n> Tests	 Manufacturer Fabricant Hersteller
	Consult instructions for use Voir les instructions d'utilisation Gebrauchsanweisung beachten	LOT Batch code Numéro de lot Chargenbezeichnung	